



REPOSITIONING AIDS

SLIDE SHEETS - FOR LATERAL TRANSFER AND REPOSITIONING

Reusable



- Washable
- Low-friction nylon
- Tubular design

Disposable



50 Sheet dispenser pack, Part # 2.01
(31" x 13" x 5"; 13 lbs.
Each sheet is 78" x 30")
Optional wall mount bracket Part # 2.09
Optional door mount bracket Part # 2.10

- No laundry required
- Wall and door mount options
- Non-woven fabric
- Low-friction on one side
- Instructions printed on dispenser pack
- Single patient use for maximum infection control

GLIDE BOARDS - FOR SEATED TRANSFER

Self-Glide
Part # 470.0050
20" long

Tetra Glide
Part # 470.0863
25" long


Curved
Part # 470.0012
27" long

Great
Part # 470.0080
31" long

Biggest
Part # 470.0102
40" long



Weight limit: 300 lbs.



Tetra Glide model has unique shape to "hook" around wheel chair arm.

- Ergonomically tapered ends
- Rubber grips on bottom
- Low-friction, medical grade plastic
- Easy to clean
- Built-in handle



TURN PLATE - FOR STANDING PIVOT



Part # 470.0380
15" Diameter



- Easily hangs on wall or door • Medical grade plastic
- Rubber grip on bottom • Use on any floor surface
- Easy to wipe clean

ROTATING SEAT CUSHION



Large
Part # 8.01
18" Diameter
(shown)

Small
Part # 460.1038
15" Diameter

- Can be used in car, chair or bed • Breathable, wickable material
- Soft and comfortable • Comes apart easy for cleaning • Can be laundered
- 2-piece construction with low-friction interior surfaces

ONE-WAY SLIDE CUSHION

- Slides back but NOT forward
- Maintains safe, secure posture
- Prevents fall-related injury
- Clearly labeled for proper orientation
- Soft, comfortable breathable/wickable material; tubular design
- Can be used in car, chair or bed



Part # 450.4060
18" x 16"

